

Laboratory Explosion Proof & Re-Circulating Cooling Freezer

실험실용 내부 방폭형 냉동고

배송비 별도

Cat. No.	Model	용량	온도 범위	EA Price ₩
F06-540-015	FROSTER-LABEX-330	300 L	approx. -5 ~ -30°C	18,155,100
F06-540-020	FROSTER-LABEX-530	500 L		19,217,900
F06-540-025	FROSTER-LABEX-730	700 L		20,856,500

- 전기 안전 인증 KC 인증 획득
- 안정제어
 - 안정성을 최대화 하기 위해 Micro Processor에 의해 온도 제어
 - 전원 장애 또는 오작동의 경우 시각/청각 경보기 울림
- 견고한 하우징 (아연 도금 강판 및 녹방지 코팅)
- 내부는 불연성 재질로 되어 있음
- 선반 높이 자유롭게 조절 가능
- 강제공기 냉각 시스템으로 온도를 제어 및 균일한 온도 분포 유지
- 빠른 자동 제상
- 저소음 Compressor
- EMC Directive 2004/108/EC (유럽 CE 전자파적합성기기 인증)
- Validation Service (IQ, OQ, PQ) 가능
- Kirsch

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
F06-540-540	SHB-FEX330	Shelf including Bracket for LABEX-330	202,300
F06-540-544	SHB-FEX530	Shelf including Bracket for LABEX-530	202,300
F06-540-548	SHB-FEX730	Shelf including Bracket for LABEX-730	202,300
F06-540-555	SHR-FEX330	Shelf including Rail for LABEX-330	357,000
F06-540-559	SHR-FEX530	Shelf including Rail for LABEX-530	357,000
F06-540-563	SHR-FEX730	Shelf including Rail for LABEX-730	357,000
F06-540-570	LD-FEX530	Length divider for LABEX-530	26,300
F06-540-575	CD-FEX530	Cross divider for LABEX-530	8,600
F06-540-580	WB-FEX330	Wire Basket including Rail for LABEX-330	555,100
F06-540-584	WB-FEX530	Wire Basket including Rail for LABEX-530	555,100
F06-540-588	WB-FEX730	Wire Basket including Rail for LABEX-730	555,100



FROSTER-LABEX-330



FROSTER-LABEX-530



FROSTER-LABEX-730

Cat. No.	F06-540-015	F06-540-020	F06-540-025	
Model	FROSTER-LABEX-330	FROSTER-LABEX-530	FROSTER-LABEX-730	
Capacity	300 L	500 L	700 L	
Temperature Setting	approx. -5 ~ -30°C (±2.0°C)			
Voltage	220~240V, 50Hz			
Power Consumption	577 W	510 W		
Normal Consumption (kWh/24h)	3.98	4.70	5.70	
Admissible Ambient Temperature	10 ~ 32°C			
Heat Emission (Max.)	711 W	690 W		
Material	External	Galvanised Sheet Steel with White, Anti-Scratch Powder Coating		
	Internal	Smooth Aluminium with Colorless Protective Coating		
	Shelf	Plastic Coated		
Dimension (cm)	Overall	w74×d78×h159	w77×d76×h193~196	w77×d98×h193~196
	Interior	w59×d52×h95	w62×d57×h140	w62×d77×h140
	Exterior (w door open at 90°)	w74×d142×h159	w77×d144×h193~196	w77×d166×h193~196
	Shelf	w57×d42	w59×d45	w59×d65
개수	선반	4	5	
Max. Load (kg)	선반	40		
무게 (kg)	Net	120	140	165
	Gross	135	170	200