

Vacuum Cold Trap Bath 진공 저온 동결 트랩

배송비 별도

Cat. No.	Model	용량	Temp.	EA Price ₩
B04-01-310	VCT42	4.2 L	Down to -40°C	3,604,000
B04-01-320	VCT85		Down to -70°C	6,107,200

- 진공 라인에서 발생하는 증기 및 유해성 물질 응축용
- 열 매체유를 사용하지 않고, Cooling Coil을 이용한 직냉방식으로 용매를 액화 시킴
- Trap 되어지는 양이 많을 때 사용이 편리함
- Glass Trap을 사용하지 않아 안전하고 편리함
- 온도 조절 불가능
- 컴팩트 저소음형 동결 트랩
- 착탈식 전면 커버 구조, 냉동기 유지 관리 용이, 효율 증대
- Bath 내부 용액을 편리하게 배출 시킬 수 있는 Drain 장착
- Graphic LCD Display, 기능간편하고 직관적인 Controller 조작
- 온도 Calibration 기능
- Automatical Fault Sensing System 기기 스스로 이상 상황을 감시 (특허 등록)
- **PL Insurance (제조물 배상 책임 보험)**
- LK Labkorea

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
B04-01-A40	GCT02	Spare Rubber Gasket	65,800



Drain



Bath Cover



Rubber Gasket

Cat. No.	B04-01-310		B04-01-320	
Model	VCT42		VCT85	
Capacity	4.2 L			
Controller	Control	Digital Control		
	Resolution	0.1°C		
Temperature	Range	Down to -40°C	Run start / Wait start (99 hr 59 min)	
	Accuracy	± 0.3°C		
Size	Interior	Ø210×h220 mm		
	Exterior	w340×d460×h640 mm	w450×d610×h750 mm	
Material	Interior	Stainless Steel 304		
	Exterior	Powder Coated Steel		
	Cover	Stainless Steel 304		
Electric Supply	Gasket	Sillicon Rubber		
	Power	1 Phase, 220 Vac, 60 Hz		
Other	Max Consumption	400 W (1.9 A)	800 W (3.7 A)	
	Power Line	Standard Plug		
Safety	Drain	1/4" Hose Nipple Male		
	Caster	Leveling Caster		
Safety		Over Current Protection, Automatical Fault Sensing System		

Vacuum Oven은 744 ~ 748 Page 참조

