

## Modular HDPE Safety Cabinet HDPE 조립식 안전 캐비닛, 산 및 부식성용

배송비 별도

Cat. No.	Model	용량	선반 수	Size, wxdxh	EA Price ₩
S01-146-170	M04CRA	15 L	2	457×457×559mm	777,600
S01-146-174	M12CRA	45 L		457×559×914mm	855,400
S01-146-178	M24CRA	90 L	4	M12CRA 2개	1,710,700
S01-146-182	M48CRA	181 L	8	M12CRA 4개	3,110,400

- Modular Quick-Assembly® Poly Cabinet
- 100% HDPE 구조로 질산, 황산, 염산과 같은 강산 저장용으로 알맞음
- 조립 도구 없이 간단하게 조립이 가능함
- M04CRA : 4 Gallon (15 L) Stand Alone Unit
- M12CRA : 12 Gallon (45 L) 제품으로 추가하여 M24CRA나 M48CRA로 조립 가능
- Lockable Door는 180도 오픈 가능
- 두 개의 2" Poly Vent와 Cap
- Poly 선반의 높이 조절 가능 (3단계)
- 모듈 조립을 위한 Bracket은 포함되어져 있음
- Eagle



S01-146-174



S01-146-174 Open 모습



수직 결합 S01-146-178



수평 결합 S01-146-178

## Accessories for "Eagle" Safety Cabinets EAGLE 안전 캐비닛 액세서리

- Eagle

### Metal Shelf for One Door Safety Cabinet

Cat. No.	Model	Description	EA Price ₩
S01-146-918	29935	15L 용 (1904X, 1903X, CRA-1904X, CRA-1903X)	105,300
S01-146-919	29936	45L, 60L, 90L 용 (1925X, 1906X 1923X, 1924X, 1925X, 2310X, CRA-1923X, CRA-2310X)	125,900
S01-146-923	29941	181L 용 (1946X)	240,400



29935

### Metal Shelf for Two Door Safety Cabinet

Cat. No.	Model	Description	EA Price ₩
S01-146-920	29937	56L, 113L, 170L 용 (ADD-15X, ADD-14X, 1932X, 1947X, 1930X, 1945X, 3010X, 4510X, CRA-32X, CRA-47X, CRA-3010X, CRA-4510X, CRA-30X, CRA-45X)	175,500
S01-146-921	29938	90 L 용 (1976X, 1975X)	145,100
S01-146-922	29939	83 L 용 (1971X, 1970X, CRA-71X, CRA-70X)	221,300
S01-146-924	29944	227 L 용 (1962X, 6010X, CRA-62X, CRA-610X)	291,900
S01-146-925	29945	341 L 용 (1992X, 9010X)	396,900



29938

### Antistatic Wire

Cat. No.	Model	Description	EA Price ₩
S01-146-940	08500	Antistatic Wire with Dual Hand Clamps	75,200
S01-146-165	1910G	Self-Closing Adaptor	694,800



08500



1910G

### Spare Poly Tray for Acid & Corrosive Cabinet

Cat. No.	Model	Description	EA Price ₩
S01-146-960	29050	15L 용 (CRA-1903X, CRA-1904X)	184,800
S01-146-962	29051	90L 용 (CRA-1923X, CRA-2310X)	362,400
S01-146-964	29052	83L 용 (CRA-71X, CRA-70X)	146,100
S01-146-966	29053	113, 170 L 용 (CRA-32X, CRA-47X, CRA-4510X, CRA-30X, CRA-45X)	260,000
S01-146-968	29054	227 L 용 (CRA-62X, CRA-6010X)	344,000



29050, 29051



29053