

Deck Mounted Eyewash & Drench Hose 고정형 세안기 & 드렌치 호스

Cat. No.	Model	Type	Hose	EA Price ₩
S01-518-450	G5022	Deck Mounted	8피트 PVC 강화 호스 최대 압력 : 275 PSI	1,405,100

- 싱크대, 실험대 판에 고정하여 눈세척
- 드렌치 호스로 눈, 얼굴, 신체 세척
- 스테인레스 스틸 레버 핸들에 의해 밸브 작동됨
- 스퀴즈 핸들 고정클립을 사용하면, 스퀴즈 핸들을 잡지 않고도 물이 분사됨
- 내부에 듀얼 역류방지 밸브 내장
- ANSI / ISEA Z358.1 - 2014 준수

ANSI(American National Standards Institute) : 미국 표준 협회
ISEA(International Safety Equipment Association) : 국제 안전 장비 협회

- Guardian



Safety Station with Eyewash 세안기 결합형 비상 샤워기

- 비상 샤워기에 세안기(눈 세척 용도)가 결합된 형태
- 샤워 헤드 재질 : 10" 오렌지색 ABS 플라스틱
- G1902BC는 먼지나 이물질 등으로부터 오염되지 않도록 스테인레스 스틸 커버 장착
- G1902HFC는 작동판(Flag Handle)이나 발판(Foot treadle)으로 작동
- G1990, G1992는 강산, 부식성 물질을 취급하는 장소에 적합
- G1941은 동파방지 밸브 장착 (물 온도가 4°C 이하 일 때 밸브가 자동으로 열려 물이 순환하여 물이 얼지 않음)
- ANSI / ISEA Z358.1 - 2014 준수
- ANSI(American National Standards Institute) : 미국 표준 협회
- ISEA(International Safety Equipment Association) : 국제 안전 장비 협회
- Guardian

기본형

Cat. No.	Model	Type	Bowl 재질	EA Price ₩
S01-518-337	G1902P	Safety Station	ABS Plastic	2,924,300
S01-518-340	G1902	Safety Station with Eyewash	Stainless Steel	3,260,800

스테인레스 스틸 보울&커버

Cat. No.	Model	Type	Cover 재질	EA Price ₩
S01-518-343	G1902BC	Safety Station with Eyewash, w/ Cover	Stainless Steel	4,069,400

작동판&발판 작동형

Cat. No.	Model	Type	Bowl 재질	EA Price ₩
S01-518-346	G1902HFC	Safety Station with Eyewash, w/ Flag handle & Foot treadle	Stainless Steel	3,517,300

PVC 재질

Cat. No.	Model	Type	재질	EA Price ₩
S01-518-349	G1990	Safety Station with Eyewash	PVC, PVC Coating (부식 방지), Stainless Steel ball valve	5,946,500
S01-518-352	G1992	Safety Station with Eyewash	PVC, PVC Coating (부식 방지)	5,461,000

동파방지밸브 장착형

Cat. No.	Model	Type	Bowl 재질	EA Price ₩
S01-518-355	G1941	Safety Station with Eyewash, w/ Freeze Protection Valve	Stainless Steel	6,197,000



G1902P



G1902



G1902HFC




G1992

Emergency Shower & Eye Wash 비상 샤워기와 세안기 일체형

Cat. No.	Model	Type	Inlet Pipe	EA Price ₩
S01-121-017	SES-0010	Red Epoxy ③	38 mm	893,300
S01-121-019	SES-0013	Polished ④	51 mm	956,600

- 재질 : Stainless Steel 304
- Emergency Shower : Ø200 mm
- Eye Wash Bowl : 300 mm
- Eye Wash Spray Head : PP, Built-in Stainless Filter
- LabSafety®



③



④