

Safety Waste System, Explosion Protected, SCAT® 안전 폐기물 처리 시스템

Cat. No.	Model	용량	Funnel top	Thread	EA Price ₩
B08-1034-100	307320	10 L	ID 200 mm	S60/61	989,500
B08-1034-104	307321	20 L			1,001,000
B08-1034-110	307416	10 L with Tray			1,784,500
B08-1034-114	307315	20 L with Tray			1,307,500
B08-1034-120	117985	Electrostatic conductive Tray HDPE			1,069,200

- 재질 : Electroconductive HDPE (High Density Polyethylene)
- Intergrated Level Control 포함
- 상부 Funnel에는 탈착 가능한 Sieve가 부착되어 있어 찌꺼기를 걸러줌 또한 탈착이 가능하여 세척에도 용이함
- 폐기용액의 임시 보관용 (Etc, Solvent Waste...)
- 구성
 - 1) Safety Electroconductive HDPE Container
 - 2) Safety Electroconductive HDPE Funnel (Ø200 mm)
 - 3) Antistatic Wire
- SCAT Europe

Safety HDPE Container Only

Cat. No.	Model	용량	Thread	EA Price ₩
S01-1034-006	108042	10 L	S60/61	433,300
S01-1034-012	108043	20 L		552,400

- 재질 : Electroconductive HDPE (High Density Polyethylene)
- Floating Level Control을 통하여 용기 내 폐액량 확인 가능
- Electrically Conductive 제품의 Antistatic Wire는 별매

Safety Waste Funnel Only

Cat. No.	Model	Funnel Top	Thread	EA Price ₩
S01-1034-406	318980	140 mm	S60/61	231,400
S01-633-470	317621	200 mm		465,200

- 재질 : Electroconductive HDPE (High Density Polyethylene)
- Antistatic Wire & Removable Sieve 포함
- Funnel은 Self-Closing Ball Valve가 장착되어 있음
- 하부 Stem 길이 : 220 mm

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
S01-1034-505	318989	Funnel용 Spare Sieve, Ø140 mm	28,700
S01-633-821	108009	Grounding cable, 2 clamps, L = 1,5 m	109,400



307416



307321



Grounding Cable

Foldable Canister / Politainer, SCAT® 접이식 [공간절약형] 폐액통

Cat. No.	Model	용량	Size (w×d×h)	Thread	EA Price ₩
B08-1034-010	107330	5 L	178×178×178 mm	GL38	22,000
B08-1034-020	107331	10 L	228×228×228 mm		26,500
B08-1034-030	107333	20 L	285×285×285 mm		S60/61

- 비사용시 접어서 보관이 가능하여 공간 절약할 수 있음
- 편리한 사후처리 가능
 - UN-Y 위험물 처리 인증 (Group II & III)
- SCAT Europe

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
S01-1034-910	107334	Box for Politaner 5 L	17,600
S01-1034-920	107335	Box for Politaner 10 L	35,400
S01-1034-930	107336	Box for Politaner 20 L	39,300
S01-1034-990	107338	PTFE Thread Adapter, GL38 → GL45	200,000
S01-1034-995	107339	PP Thread Adapter, GL38 → GL45	111,100

- Thread Adapter 사용시 상부 Funnel 장착하여 사용 가능함

