

I Accessory for IKA Circulator
IKA 항온 순환기 약세사리

For Immersion Circulator

Plastic Water Bath, Max. Temp 100°C

Cat. No.	Model	Capacity	Size	EA Price ₩
B04-05-111	4248100	8 L	w320×d335×h155 mm, S	510,000
B04-05-113	4248200	18 L	w338×d584×h155 mm, L	470,000
B04-05-115	20006884	24 L	w338×d584×h205 mm, L	1,350,000
B04-05-123	20006341	25 L	w365×d722×h165 mm, L	2,190,000
B04-05-125	20004382	9 L	w287×d356×h167 mm, S	690,000
B04-05-127	20004383	20 L	w365×d553×h167 mm, L	920,000
B04-05-129	20007044	52 L	w365×d648×h316 mm L	1,260,000

Stainless Steel Bath, Max. Temp 200°C

Cat. No.	Model	Capacity	Size	EA Price ₩
B04-05-117	4248500	9 L	w374×d377×h195 mm, S	1,070,000
B04-05-119	4577500	12 L	w354×d461×h195 mm, M	1,270,000
B04-05-121	4248600	20 L	w354×d641×h195 mm, L	1,410,000

Bath Cover

Cat. No.	Model	Description	적용 모델	EA Price ₩
B04-05-141	4471500	Small Cover, S	ICC	250,000
B04-05-143	25000290	Medium Cover, M		310,000
B04-05-145	4471600	Large Cover, L		410,000
B04-05-131	4577600	Medium Cover, L	IC	210,000
B04-05-133	4471800	Large Cover, L		250,000

Bath Bridge

Cat. No.	Model	Description	적용 모델	EA Price ₩
B04-05-147	20003077	Small Bridge, S	ICC	190,000
B04-05-135	4472800	Small Bridge, M	IC	210,000
B04-05-149	20003078	Large Bridge, L	ICC	220,000

Cooling Coil and Pump Connection Set


Cat. No.	Model	Description	적용 모델	EA Price ₩
B04-05-137	20005116	CC1 Cooling Coil	ICC	130,000
B04-05-151	25001061	CC2 Cooling Coil	IC	190,000
B04-05-153	4471900	PCS.ICC Pump connection set	ICC	260,000



Plastic Bath



Stainless Steel Bath



Bath Cover



Bath Bridge



CC1 Cooling Coil



CC2 Cooling Coil



Pump Connection Set



Expandable with Inlay 1~3



Variable Rack



Floating Tube Rack

Bath Insert

Cat. No.	Model	Description	Tube Size	Immersion Depth	Bath Size	적용 모델	EA Price ₩
B04-05-185	20004026	Stainless Steel Tube Rack, S	Ø13 mm×57 hole	70 mm	Small	ICC, IC 공용	210,000
B04-05-187	20004027		Ø17 mm×37 hole	100 mm			210,000
B04-05-189	20004028		Ø22 mm×22 hole	50 mm			210,000
B04-05-191	20004029	Stainless Steel Tube Rack, M/L	Ø13 mm×73 hole	70 mm	Medium, Large		360,000
B04-05-193	20004030		Ø17 mm×47 hole	100 mm			360,000
B04-05-195	20004031		Ø22 mm×30 hole	50 mm			360,000
B04-05-197	20004032	Expandable with Inlay 1 ~ 3			0 ~ 120 mm		460,000
B04-05-199	20004033	Variable Rack, Inlay 1, 2 pcs	Ø13 mm×84 hole	All Size			150,000
B04-05-201	20004034	Variable Rack, Inlay 2, 2 pcs	Ø17 mm×51 hole				150,000
B04-05-203	20004035	Variable Rack, Inlay 3, 2 pcs	Ø22 mm×33 hole		150,000		
B04-05-211	20003667	Floating Tube Rack, 5 pcs	1.5/2.0 ml×24 hole		All Size	85,000	
B04-05-213	20003668		15 ml×8 hole			85,000	
B04-05-215	20003669		50 ml×4 hole			85,000	
B04-05-205	20004614	Stainless Steel Variable Bottom for size L				Large	910,000
B04-05-207	20006212	Bottom Plate for Bath size L					1,090,000
B04-05-209	20007353	Bottom Plate for Bath size L					910,000

다음 페이지에서 계속...