

## IKA Small Orbital Shaker 소형 오비탈 셰이커, KS 130 basic & control

Cat. No.	Model	Type	Display	EA Price ₩
S04-05-107	2980000	KS 130 basic	LED Line	2,670,000
S04-05-114	2980100	KS 130 control	Digital	4,610,000

- "Circular 원형운동"의 소형 Shaking & Vibration system, Culture & Mixing용
- 무소음, 무부하의 균일 운동, 고장이 없는 반영구적인 수명
- 다용도, 다양성 : 25 ~ 500 ml의 어떤 모양이나 크기의 vessel도 수용 가능
- 기본 Unit에 Platform등 각종 부품은 별도
- KS130 Control은 PC와 연결 상태 저장 및 작동 할 수 있음
- Protection class acc. to DIN EN 60529 (IP21)
- IKA

Cat. No.	S04-05-107	S04-05-114
Type	KS 130 basic	KS 130 control
Shaking Movement	Orbital	
Orbital Diameter	4 mm	
Motor Rating Input	45 W	
Motor Rating Output	10 W	
Speed Range	0 ~ 800 rpm	
Timer	∞ / 5 ~ 50 min	∞ / 0 ~ 9 hr 59 min
Size	w270×d316×h98 mm	
Interface	-	RS 232
Power	230V, 50/60 Hz	

### Optional Accessory


Cat. No.	Model	Type	Description	EA Price ₩
S04-05-506	8017300	AS130.1	Universal Attachment	1,070,000
S04-05-512	3115000	AS130.2	Fixing Clip Attachment	340,000
S04-05-518	3120000	AS130.3	Dish Attachment	260,000
S04-05-514	3120300	AS130.4	Test Tube Attachment	600,000
S04-05-912	1234300	AS2.1	Fixing Clip, 25 ml용	60,000
S04-05-918	1234400	AS2.2	Fixing Clip, 50 ml용	60,000
S04-05-924	1234500	AS2.3	Fixing Clip, 100 ml용	60,000
S04-05-930	1234600	AS2.4	Fixing Clip, 200 / 250 ml용	120,000
S04-05-936	1234700	AS2.5	Fixing Clip, 500 ml용	120,000
H03-05-650	4720000	labworldsoft	Software	10,190,000

- AS 130.1 Universal Attachment : w325 × d234 × h88 mm
  - For use with various types of vessels
  - Set-up Plate : w220 × d220 mm
- AS 130.2 Fixing Clip Attachment : w230 × d230 × h24 mm
  - For use with round flasks, erlenmeyer flasks
  - AS2.1 : 25 ml용, 20개 장착 가능
  - AS2.2 : 50 ml용, 12개 장착 가능
  - AS2.3 : 100 ml용, 12개 장착 가능
  - AS2.4 : 250 ml용, 4개 장착 가능
  - AS2.5 : 500 ml용, 4개 장착 가능
- AS 130.3 Dish Attachment : w420 × d270 × h32 mm
  - For Smooth Shaking Operations in the Low Viscosity
  - For Use Petri Dishes or Culture Bottles
  - PP Slip-resistant Foil 포함
  - Set-up Plate : w220 × d340 mm
- AS 130.4 Test tube Attachment : w220 × d230 × h95 mm
  - For use with Test Tube, Small Tube, Centrifuge Tube
  - Ø10~16 mm, 64 Hole
  - 최소 Vessel 높이 : 80 mm



CE

KS 130 basic



KS 130 control




AS 130.1



AS 130.4



AS 130.3



AS 130.2



Fixing Clip