

Accessory for IKA Circulator 에서 계속...

## For All IKA Circulator

### Hose

Cat. No.	Model	Description	EA Price ₩
B04-05-165	4568800	PVC Tubing 2 pcs id 8mm, od 12mm, 길이 1.5 m	35,000
B04-05-167	4568900	PVC Tubing 2 pcs id 12mm, od 16mm, 길이 1.5 m	45,000
B04-05-169	4569000	Silicone Tubing 2 pcs id 8mm, od 12mm, 길이 1.5 m	85,000
B04-05-171	4569100	Silicone Tubing 2 pcs id 12mm, od 16mm, 길이 1.5 m	100,000
B04-05-177	4569200	Viton(FKM), Stainless Steel Tubing 2 pcs id 8mm, od 12mm, 길이 1.5 m	730,000
B04-05-179	4569300	Viton(FKM), Stainless Steel Tubing 2 pcs id 10mm, od 16mm, 길이 1.5 m	960,000
B04-05-173	20004612	PUR Tubing 2 pcs id 8mm, od 12mm, 길이 1.5 m	150,000
B04-05-175	20004613	PUR Tubing 2 pcs id 12mm, od 16mm, 길이 1.5 m	180,000
B04-05-181	4569400	Insulation hose 2 pcs id 13mm, thickness 8mm	70,000
B04-05-183	4569500	Insulation hose 2 pcs id 17mm, thickness 12mm	90,000
B04-05-168	20102324	Silicone Tubing id 6mm, od 9mm, 길이 3 m	65,000
R02-05-311	LT 5.30	Tempering hose 800 mm, 2pcs	1,330,000
R02-05-314	LT 5.31	Tempering hose 1200 mm, 2pcs	1,480,000
R02-05-317	LT 5.32	Tempering hose 1600 mm, 2pcs	1,630,000

### Optional Valve

Cat. No.	Model	Description	EA Price ₩
B04-05-217	20003763	MV 1 Magnetic Valve	990,000
B04-05-219	20000249	CO V 1 Closed pressure Valve	1,450,000
B04-05-221	20004620	Ball Valve M16×1	210,000
B04-05-223	20004618	Mechanic fluid level controller	750,000

### Optional Adapter

Cat. No.	Model	Description	EA Price ₩
B04-05-225	25001212	Elbow connector 90°, M16×1 male thread	200,000
B04-05-227	20004667	id 6mm용 Olive connector, 2 pcs	70,000
B04-05-229	20004566	id 8mm용 Olive connector, 2 pcs	70,000
B04-05-231	20004568	id 10mm용 Olive connector, 2 pcs	70,000
B04-05-233	20004889	id 12mm용 Olive connector, 2 pcs	70,000
B04-05-235	20004569	NPT 1/4, 2 pcs	110,000
B04-05-237	20004570	NPT 1/2, 2 pcs	310,000
B04-05-239	20004571	NPT 3/4, 2 pcs	390,000
B04-05-241	20004583	Lock Nut, 2 pcs	85,000
B04-05-243	20004584	Stopper, 2 pcs	85,000

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
H03-05-651	20017366	Labworldsoft® 6 pro	23,240,000
H03-05-650	4720000	Labworldsoft® 6 Starter	10,190,000
H03-05-664	20105873	Labworldsoft® 6 Advanced	17,750,000
B04-05-245	20000984	WH10 WiCo Holder	55,000
B04-05-247	4284700	PT 100.30 Temp. Sensor, 250 mm	740,000
B04-05-249	20004629	PT 100 extension (LEMO)	460,000
B04-05-251	20003666	PP Floating Globe, Ø20 mm, 500 pcs	280,000
B04-05-155	20003521	Silicone Oil, 10 kg, Red	1,240,000
B04-05-157	20003520	Silicone Oil, 10 kg, Clear	1,120,000
B04-05-159	20003518	Silicone Oil, 9 kg	1,170,000
C30-05-718	20014031	Water protect, 100 ml	180,000



B

Bath

Accessory for IKA Circulator