

IKA Orbital Shaker 중형 오비탈 셰이커, KS 260 basic & control

Cat. No.	Model	Type	Display	EA Price ₩
S04-05-121	2980200	KS 260 basic	LED Line	3,600,000
S04-05-128	2980300	KS 260 control	Digital	6,270,000

- "Circular 원형운동"의 중형 Shaking & Vibration system, Culture & Mixing용
- 무소음, 무부하의 균일 운동, 고장이 없는 반영구적인 수명
- 다용도, 다양성 : 25 ~ 500 ml의 어떤 모양이나 크기의 vessel도 수용 가능
- 기본 Unit에 Platform 등 각종 부품은 별도
- KS206 Control은 PC와 연결 상태 저장 및 작동 할 수 있음
- Protection class acc. to DIN EN 60529 (IP21)
- IKA

Cat. No.	S04-05-121	S04-05-128
Type	KS 260 basic	KS 260 control
Shaking Movement	Orbital	
Orbital Diameter	10 mm	
Motor Rating Input	45 W	
Motor Rating Output	10 W	
Speed Range	20 ~ 500 rpm	10 ~ 500 rpm
Timer	∞ / 5 ~ 50 min	∞ / 0 ~ 9 hr 59 min
Size	w360 × d420 × h98 mm	
Interface	-	RS 232
Power	230V, 50/60 Hz	

Optional Accessory


Cat. No.	Model	Type	Description	EA Price ₩
S04-05-521	8017400	AS260.1	Universal Attachment	1,710,000
S04-05-528	3115500	AS260.2	Fixing Clip Attachment	410,000
S04-05-535	3120600	AS260.3	Dish Attachment	210,000
S04-05-542	3120900	AS260.5	Separating Funnel Attachment	1,340,000
S04-05-912	1234300	AS2.1	Fixing Clip, 25 ml용	60,000
S04-05-918	1234400	AS2.2	Fixing Clip, 50 ml용	60,000
S04-05-924	1234500	AS2.3	Fixing Clip, 100 ml용	60,000
S04-05-930	1234600	AS2.4	Fixing Clip, 250 ml용	120,000
S04-05-936	1234700	AS2.5	Fixing Clip, 500 ml용	120,000
S04-05-942	3819300	AS2.6	Fixing Clip, 1 L용	120,000
H03-05-650	4720000	labworldsoft	Software	10,190,000

- AS 260.1 Universal Attachment : w425 × d335 × h135 mm
 - For use with various types of vessels
 - Set-up Plate : w320 × d320 mm
- AS 260.2 Fixing Clip Attachment : w330 × d330 × h24 mm
 - For use with round flasks, erlenmeyer flasks
 - AS2.1 : 25 ml용, 56개 장착 가능
 - AS2.2 : 50 ml용, 23개 장착 가능
 - AS2.3 : 100 ml용, 23개 장착 가능
 - AS2.4 : 250 ml용, 11개 장착 가능
 - AS2.5 : 500 ml용, 9개 장착 가능
 - AS2.6 : 1 L용, 5개 장착 가능
- AS 260.3 Dish Attachment : w410 × d370 × h32 mm
 - For Smooth Shaking Operations in the Low Viscosity
 - For Use Petri Dishes or Culture Bottles
 - PP Slip-resistant Foil 포함
 - Set-up Plate : w320 × d320 mm
- AS 260.5 Separating Funnel Attachment : w334 × d425 × h145 mm
 - For use with Separating Funnel
 - 50 ml용, 6개 장착 가능
 - 100 ml용, 5개 장착 가능
 - 250 ml용, 3개 장착 가능
 - 500 ml용, 3개 장착 가능



CE

KS 260 basic



CE

KS 260 control




AS 260.1



AS 260.5



AS 260.3



AS 260.2



Fixing Clip