

## IKA Digital Orbital & Horizontal Shaker 대형 오비탈 & 호리존탈 셰이커

Cat. No.	Model	Type	운동	EA Price ₩
S04-05-135	2526400	KS 501 digital	Orbital	3,025,000
S04-05-142	2527000	HS 501 digital	Horizontal	3,020,000

- Orbital (원형운동), Horizontal (좌우 왕복운동)
- 대형 Shaking, Vibration system, Culture & Mixing용
- 무소음, 무부하의 균일 운동, 고장이 없는 반영구적인 수명
- 다용도, 다양성 : 25 ~ 500 ml의 어떤 모양이나 크기의 vessel도 수용 가능
- 기본 Unit에 Platform등 각종 부품은 별도
- Protection class acc. to DIN EN 60529 (IP21)
- IKA

Cat. No.	S04-05-135	S04-05-142
Type	KS 501 digital	HS 501 digital
Shaking Movement	Orbital	Reciprocating
Shaking Diameter	30 mm	
Motor Rating Input	70 W	
Motor Rating Output	19 W	
Speed Range	0 ~ 300 rpm	
Timer	∞ / 1 ~ 56 min	
Size	w505×d585×h120 mm	
Interface	RS 232	
Power	230V, 50/60 Hz	

### Optional Accessory

Cat. No.	Model	Type	Description	EA Price ₩
S04-05-549	8000200	AS501.1	Universal Attachment	2,290,000
S04-05-556	2341100	AS501.4	Fixing Clip Attachment	2,740,000
S04-05-563	2339600	AS501.5	Dish Attachment	590,000
S04-05-570	8000300	AS501.2	Separating Funnel Attachment	3,030,000
S04-05-577	8000400	AS501.3	Separating Funnel Attachment	2,580,000
S04-05-912	1234300	AS2.1	Fixing Clip, 25 ml용	60,000
S04-05-918	1234400	AS2.2	Fixing Clip, 50 ml용	60,000
S04-05-924	1234500	AS2.3	Fixing Clip, 100 ml용	60,000
S04-05-930	1234600	AS2.4	Fixing Clip, 250 ml용	120,000
S04-05-936	1234700	AS2.5	Fixing Clip, 500 ml용	120,000
S04-05-942	3819300	AS2.6	Fixing Clip, 1 L용	120,000
H03-05-650	4720000	labworldsoft	Software	10,190,000

- AS 501.1 Universal Attachment : w480×d500×h120 mm
  - For use with various types of vessels
  - Set-up Plate : w420×d420 mm
- AS 501.4 Fixing Clip Attachment : w475×d460×h95 mm
  - For use with round flasks, erlenmeyer flasks
  - AS2.1 : 25 ml용, 110개 장착 가능
  - AS2.2 : 50 ml용, 55개 장착 가능
  - AS2.3 : 100 ml용, 35개 장착 가능
  - AS2.4 : 250 ml용, 16개 장착 가능
  - AS2.5 : 500 ml용, 12개 장착 가능
  - AS2.6 : 1 L용, 5개 장착 가능
- AS 501.5 Dish Attachment : w450×d450×h45mm
  - For Smooth Shaking Operations
  - For Use Petri Dishes or Culture Bottles
  - Set-up Plate : w420×d420 mm
- AS 501.2 Separating Funnel Attachment : w480×d505×h190 mm
  - For use with Separating Funnel
  - 50 ml용, 12개 장착 가능
  - 100 ml용, 10개 장착 가능
  - 250 ml용, 6개 장착 가능
- AS 501.3 Separating Funnel Attachment
  - For use with Separating Funnel
  - 500 ml용, 4개 장착 가능
  - 1000 ml용, 3개 장착 가능
  - 2000 ml용, 2개 장착 가능

