

IKA Rocking Shaker 2D 락커, Rocker 2D Series

Cat. No.	Model	Type	Angle	EA Price ₩
S04-05-187	4002000	Rocker 2D basic	8°	1,430,000
S04-05-194	4003000	Rocker 2D digital	0 ~ 15°	1,690,000

- Touch Keypad를 이용하여 쉽게 조작 가능 • 시소 궤적으로 Rocking
- Dual Platform Option을 통하여 2가지 약세사리를 동시 사용 가능
- Speed 조절 가능 • Protection Class Acc. to DIN EN 60529: IP21
- IKA

Cat. No.	S04-05-187	S04-05-194
Type	Rocker 2D basic	Rocker 2D digital
Max. Shaking weight	2 kg	
Motor 입출력	16/9 W	
Speed Range	5 ~ 80 rpm	
Speed Display	N/A	가능
Timer	N/A	1 sec ~ 99 hr 55 min
Dimensions	w280 × d330 × 150mm	
Weight	2.2 kg	
전기 규격	110 ~ 240 V, 50/60 Hz	
DC Voltage	24 V	

Optional Accessory

Cat. No.	Model	Type	Application	EA Price ₩
S04-05-770	4532000	RD1 Mat	Universal	30,000
S04-05-777	4532100	RD2 burled Mat	Universal	35,000
S04-05-780	4581700	A 2	Second Level Attachment	40,000
S04-05-784	4532200	TB1 Tray	5 ml(∅12 mm) × 28	25,000
S04-05-791	4532300	TB2 Tray	15 ml(∅16 mm) × 16	25,000
S04-05-798	4532400	TB3 Tray	30 ml(∅25 mm) × 14	25,000
S04-05-805	4532500	TB4 Tray	50 ml(∅28 mm) × 10	40,000

Digital Horizontal Shaker 디지털 좌우 왕복형 셰이커

Cat. No.	Model	Type	Angle	EA Price ₩
S04-05-201	4000000	Rocker 3D basic	8°	1,360,000
S04-05-208	4001000	Rocker 3D digital	0 ~ 15°	1,670,000

- Touch Keypad를 이용하여 쉽게 조작 가능
- 시소와 텀블링(복합) 궤적으로 Rocking
- Dual Platform Option을 통하여 2가지 약세사리를 동시 사용 가능
- Speed 조절 가능 • Protection Class Acc. to DIN EN 60529: IP21
- IKA

Cat. No.	S04-05-201	S04-05-208
Type	Rocker 3D basic	Rocker 3D digital
Max. Shaking weight	2 kg	
Motor 입출력	16/9 W	
Speed Range	30 rpm (Fixed)	5 ~ 80 rpm
Speed Display	N/A	가능
Timer	N/A	1 sec ~ 99 hr 55 min
Size	w280 × d330 × 165 mm	w280 × d330 × 185 mm
Weight	2.2 kg	
전기 규격	110 ~ 240 V, 50/60 Hz	
DC Voltage	24 V	

Optional Accessory

Cat. No.	Model	Type	Application	EA Price ₩
S04-05-770	4532000	RD1 Mat	Universal	30,000
S04-05-777	4532100	RD2 burled Mat	Universal	35,000
S04-05-780	4581700	A 2	Second Level Attachment	40,000
S04-05-784	4532200	TB1 Tray	5 ml(∅12 mm) × 28	25,000
S04-05-791	4532300	TB2 Tray	15 ml(∅16 mm) × 16	25,000
S04-05-798	4532400	TB3 Tray	30 ml(∅25 mm) × 14	25,000
S04-05-805	4532500	TB4 Tray	50 ml(∅28 mm) × 10	40,000

