

## IKA MATRIX Orbital Thermoshaker 에서 계속...

## Spare Accessory

## MATRIX Orbital

Cat. No.	Model	Description	EA Price ₩
H04-05-684	30001157	Insert for 0.5 ml vial	85,000
H04-05-685	30001158	Insert for 1.5/ 2.0 ml vial	85,000
H04-05-686	30001159	Insert for 0.2 ml PCR vial	85,000

## Optional Accessory

## MATRIX Orbital Delta

Cat. No.	Model	Description	EA Price ₩
H04-05-672	30001253	Cover for MATRIX Orbital	220,000
H04-05-673	30001244	Heatable Matrix ThermoCover	1,010,000

## MATRIX Orbital Delta Plus

Cat. No.	Model	Description	EA Price ₩
H04-05-674	30000853	Attachment for 96 PCR Plate	1,210,000
H04-05-675	30000857	Attachment for 384 PCR Plate	1,300,000
H04-05-676	30000840	Attachment for 24 × 0.5 ml vial	1,080,000
H04-05-677	30000829	Attachment for 24 × 1.5 ml vial	1,010,000
H04-05-678	30000834	Attachment for 24 × 2.0 ml vial	1,010,000
H04-05-679	30000859	Attachment for 8 × 5 ml vial	1,080,000
H04-05-680	30000844	Attachment for 24 × 12 mm vial	1,080,000
H04-05-681	30000821	Attachment for 8 × 15 ml vial	1,080,000
H04-05-682	30000863	Attachment for 4 × 50 ml vial	1,010,000
H04-05-683	30000848	Attachment for 24 × cryo vial	1,080,000
H04-05-687	30000849	Attachment for Microtiter Plate	1,210,000



## 96-Well Microplate High Speed Shaker

### 96 웰 플레이트용 고속 셰이커

Cat. No.	Model	Capacity	Speed	EA Price ₩
S04-480-025	LS-400M	4 Plate	200 ~ 1500 rpm	934,000

- Elisa Plate(96/384 well) 또는 24/48/96 well Tissue Culture Plate용 셰이커
- LCD display, One Touch Knob
- Brushless DC Motor, Microprocessor Controlled
- 자동 고장 감지 및 알람 기능
- LK Labkorea

Cat. No.	S04-480-025		
Model	LS-400M		
Capacity	4 Plate		
Speed Range	200 ~ 1500 rpm		
Speed Step	10 rpm		
Mixing Orbit	3 mm		
Controller	Display	LCD	
	Timer	1 min ~ 99h 59min	
Size	w280 × d270 × h110 mm		
Electric Supply	Power	100 ~ 240 V, 50/60 Hz	
	Max Consumption	45W	



- Thermo Shaker for Microplate는  
668 Page 참조

