

Standard Test Sieve

국산 시험용 망체

- 재질 : 황동 Frame, Stainless Screen Mesh
- 요청에 의해 Ø150mm, Ø200mm, Ø305mm, Ø450mm도 공급 가능함
- 요청에 의해 Stainless Frame 도 공급 가능
- LK Labkorea

Ø203 mm × h41 mm Sieve

Cat. No.	Model	ASTM 규격	ISO 규격	EA Price ₩
S19-124-004	BS6350	1/4 in	6.30 mm	81,600
S19-124-008	BS5600	3 1/2 in	5.60 mm	81,600
S19-124-012	BS5000		5.00 mm	81,600
S19-124-016	BS4750	No. 4	4.75 mm	81,600
S19-124-020	BS4000	No. 5	4.00 mm	71,300
S19-124-024	BS3350	No. 6	3.35 mm	71,300
S19-124-028	BS2800	No. 7	2.80 mm	71,300
S19-124-032	BS2360	No. 8	2.36 mm	71,300
S19-124-036	BS2000	No. 10	2.00 mm	71,300
S19-124-040	BS1700	No. 12	1.70 mm	71,300
S19-124-044	BS1400	No. 14	1.40 mm	71,300
S19-124-048	BS1180	No. 16	1.18 mm	71,300
S19-124-052	BS1000	No. 18	1.00 mm	71,300
S19-124-056	BS0850	No. 20	850 um	71,300
S19-124-060	BS0710	No. 25	710 um	67,400
S19-124-064	BS0600	No. 30	600 um	67,400
S19-124-068	BS0500	No. 35	500 um	67,400
S19-124-072	BS0425	No. 40	425 um	67,400
S19-124-076	BS0355	No. 45	355 um	67,400
S19-124-080	BS0300	No. 50	300 um	67,400
S19-124-084	BS0250	No. 60	250 um	67,400
S19-124-088	BS0212	No. 70	212 um	67,400
S19-124-092	BS0200		200 um	71,300
S19-124-096	BS0180	No. 80	180 um	71,300
S19-124-100	BS0160		160 um	71,300
S19-124-104	BS0150	No. 100	150 um	71,300
S19-124-108	BS0125	No. 120	125 um	71,300
S19-124-112	BS0106	No. 140	106 um	71,300
S19-124-116	BS0100		100 um	81,600
S19-124-120	BS0090	No. 170	90 um	81,600
S19-124-124	BS0075	No. 200	75 um	81,600
S19-124-128	BS0063	No. 230	63 um	92,700
S19-124-132	BS0053	No. 270	53 um	129,500
S19-124-136	BS0045	No. 325	45 um	141,400
S19-124-140	BS0038	No. 400	38 um	249,500
S19-124-144	BS0032	No. 450	32 um	330,300
S19-124-148	BS0025	No. 500	25 um	514,400
S19-124-152	BS0020	No. 635	20 um	1,267,600
S19-124-950	SIA001	Pan & Cover Set		45,100
S19-124-954	SIA005	중간체		110,000

굵은 골재용 시험체, Ø203 mm × h41 mm

Cat. No.	Model	ASTM 규격	ISO 규격	EA Price ₩
S19-124-304	WBS5300	2.12 in	53 mm	67,400
S19-124-308	WBS4500	1 3/4 in	45 mm	67,400
S19-124-312	WBS3750	1 1/2 in	37.5 mm	67,400
S19-124-316	WBS3150	1 1/4 in	31.5 mm	67,400
S19-124-320	WBS2650	1.06 in	26.5 mm	67,400
S19-124-324	WBS2240	7/8 in	22.4 mm	67,400
S19-124-328	WBS1900	3/4 in	19 mm	67,400
S19-124-332	WBS1600	5/8 in	16 mm	67,400
S19-124-336	WBS1320	0.53 in	13.2 mm	67,400
S19-124-340	WBS1120	7/16 in	11.2 mm	67,400
S19-124-344	WBS0095	3/8 in	9.5 mm	67,400
S19-124-348	WBS0008	5/16 in	8 mm	67,400
S19-124-352	WBS0067	0.265 in	6.7mm	67,400



Pan & Cover Set

Ø75 mm × h22 mm Sieve

Cat. No.	Model	ASTM 규격	ISO 규격	EA Price ₩
S19-124-504	BS6350-75	1/4 in	6.35 mm	68,700
S19-124-508	BS5600-75	3 1/2 in	5.60 mm	68,700
S19-124-512	BS5000-75		5.00 mm	68,700
S19-124-516	BS4750-75	No. 4	4.75 mm	68,700
S19-124-520	BS4000-75	No. 5	4.00 mm	59,600
S19-124-524	BS3350-75	No. 6	3.35 mm	59,600
S19-124-528	BS2800-75	No. 7	2.80 mm	59,600
S19-124-532	BS2360-75	No. 8	2.36 mm	59,600
S19-124-536	BS2000-75	No. 10	2.00 mm	59,600
S19-124-540	BS1700-75	No. 12	1.70 mm	59,600
S19-124-544	BS1400-75	No. 14	1.40 mm	59,600
S19-124-548	BS1180-75	No. 16	1.18 mm	59,600
S19-124-552	BS1000-75	No. 18	1.00 mm	59,600
S19-124-556	BS0850-75	No. 20	850 um	59,600
S19-124-560	BS0710-75	No. 25	710 um	57,000
S19-124-564	BS0600-75	No. 30	600 um	57,000
S19-124-568	BS0500-75	No. 35	500 um	57,000
S19-124-572	BS0425-75	No. 40	425 um	57,000
S19-124-576	BS0355-75	No. 45	355 um	57,000
S19-124-580	BS0300-75	No. 50	300 um	57,000
S19-124-584	BS0250-75	No. 60	250 um	57,000
S19-124-588	BS0212-75	No. 70	212 um	57,000
S19-124-592	BS0200-75		200 um	59,600
S19-124-596	BS0180-75	No. 80	180 um	59,600
S19-124-600	BS0160-75		160 um	59,600
S19-124-604	BS0150-75	No. 100	150 um	59,600
S19-124-608	BS0125-75	No. 120	125 um	59,600
S19-124-612	BS0106-75	No. 140	106 um	59,600
S19-124-616	BS0100-75		100 um	68,700
S19-124-620	BS0900-75	No. 170	90 um	68,700
S19-124-624	BS0075-75	No. 200	75 um	68,700
S19-124-628	BS0063-75	No. 230	63 um	78,400
S19-124-632	BS0053-75	No. 270	53 um	109,300
S19-124-636	BS0045-75	No. 325	45 um	120,000
S19-124-640	BS0038-75	No. 400	38 um	211,500
S19-124-644	BS0032-75	No. 450	32 um	280,400
S19-124-648	BS0025-75	No. 500	25 um	437,200
S19-124-652	BS0020-75	No. 635	20 um	1,076,300
S19-124-654	SIA002	Pan & Cover Set		38,000