

Stainless Steel Woven Wire Mesh

스테인레스 망


국산 Stainless Steel 망

Cat. No.	Model	Opening	Wire	Mesh No.	EA Price ₩
S19-250-102	SM-2370	2.37 mm	0.8 mm	8	62,700
S19-250-108	SM-2470	2.47 mm	0.7 mm		44,600
S19-250-114	SM-2570	2.57 mm	0.6 mm		33,800
S19-250-120	SM-2670	2.67 mm	0.5 mm		25,400
S19-250-138	SM-1840	1.84 mm	0.7 mm	10	54,300
S19-250-144	SM-1940	1.94 mm	0.6 mm		41,000
S19-250-150	SM-2040	2.04 mm	0.5 mm		30,200
S19-250-162	SM-2140	2.14 mm	0.4 mm	12	21,700
S19-250-168	SM-1510	1.51 mm	0.6 mm		48,200
S19-250-174	SM-1610	1.61 mm	0.5 mm		35,000
S19-250-180	SM-1660	1.71 mm	0.4 mm		25,400
S19-250-186	SM-1310	1.31 mm	0.5 mm	14	39,800
S19-250-192	SM-1410	1.41 mm	0.4 mm		27,800
S19-250-198	SM-1460	1.46 mm	0.35 mm		21,700
S19-250-204	SM-0980	0.98 mm	0.6 mm		63,900
S19-250-210	SM-1080	1.08 mm	0.5 mm	16	44,600
S19-250-216	SM-1180	1.18 mm	0.4 mm		31,400
S19-250-222	SM-1230	1.23 mm	0.35 mm		25,400
S19-250-228	SM-1290	1.29 mm	0.29 mm		19,300
S19-250-234	SM-0960	0.96 mm	0.5 mm	18	50,700
S19-250-240	SM-1010	1.01 mm	0.4 mm		35,000
S19-250-246	SM-1060	1.06 mm	0.35 mm		38,600
S19-250-252	SM-1120	1.12 mm	0.29 mm		20,500
S19-250-258	SM-0770	0.77 mm	0.5 mm	20	55,500
S19-250-264	SM-0870	0.87 mm	0.4 mm		37,400
S19-250-270	SM-0920	0.92 mm	0.35 mm		31,400
S19-250-276	SM-1040	1.04 mm	0.29 mm		22,900
S19-250-282	SM-0830	0.83 mm	0.23 mm	24	19,300
S19-250-288	SM-0770	0.77 mm	0.29 mm		26,600
S19-250-294	SM-0710	0.71 mm	0.35 mm		35,000
S19-250-300	SM-0617	0.617 mm	0.23 mm		24,100
S19-250-306	SM-0557	0.557 mm	0.29 mm	30	31,400
S19-250-312	SM-0564	0.564 mm	0.23 mm		25,400
S19-250-318	SM-0504	0.504 mm	0.29 mm	32	33,800
S19-250-324	SM-0455	0.455 mm	0.18 mm		21,700
S19-250-330	SM-0435	0.435 mm	0.20 mm	40	24,100
S19-250-336	SM-0385	0.385 mm	0.25 mm		33,800
S19-250-342	SM-0368	0.368 mm	0.14 mm	50	19,300
S19-250-348	SM-0328	0.328 mm	0.18 mm		26,600
S19-250-354	SM-0308	0.308 mm	0.20 mm		31,400
S19-250-360	SM-0283	0.283 mm	0.14 mm		21,700
S19-250-366	SM-0243	0.243 mm	0.18 mm	60	31,400
S19-250-372	SM-0198	0.198 mm	0.12 mm		22,900
S19-250-378	SM-0178	0.178 mm	0.14 mm	80	29,000
S19-250-384	SM-0154	0.154 mm	0.10 mm		22,900
S19-250-390	SM-0132	0.132 mm	0.08 mm	120	26,600
S19-250-396	SM-0104	0.104 mm	0.065 mm	150	36,200
S19-250-402	SM-0091	0.091 mm	0.091 mm	180	45,800
S19-250-408	SM-0077	0.077 mm	0.053 mm	200	48,200
S19-250-414	SM-0062	0.062 mm	0.040 mm	250	67,500
S19-250-420	SM-0045	0.045 mm	0.040 mm	300	72,300
S19-250-426	SM-0043	0.043 mm	0.035 mm	325	79,600
S19-250-432	SM-0034	0.034 mm	0.030 mm	400	120,500
S19-250-438	SM-00258	0.0258 mm	0.025 mm	500	361,500

• Mesh No. : 임의의 Wire 중심으로부터 1"내에 있는 체의 눈금 (Opening) 수

• Size : 1 × 1 m

• LK Labkorea



Stainless Steel Woven Wire Mesh 스테인레스

Cat. No.	Model	규격	EA Price ₩
S21-643-653	TIS-060	#4 (5.25 mm)	196,600
S21-643-657	TIS-065	#5 (4.33 mm)	136,100
S21-643-661	TIS-070	#6 (3.49 mm)	181,400
S21-643-665	TIS-075	#7 (3.06 mm)	105,900
S21-643-669	TIS-080	#8 (2.61 mm)	113,500
S21-643-673	TIS-085	#9 (2.26 mm)	105,900
S21-643-677	TIS-090	#10 (2.07 mm)	100,800
S21-643-681	TIS-095	#12 (1.65 mm)	118,600
S21-643-685	TIS-100	#14 (1.45 mm)	93,300
S21-643-689	TIS-105	#16 (1.22 mm)	108,400
S21-643-693	TIS-110	#18 (1.05 mm)	123,600
S21-643-697	TIS-115	#20 (0.98 mm)	78,100
S21-643-701	TIS-120	#24 (0.73 mm)	85,700
S21-643-705	TIS-125	#30 (0.63 mm)	78,100
S21-643-709	TIS-130	#40 (0.46 mm)	78,100
S21-643-713	TIS-135	#50 (0.33 mm)	85,700
S21-643-717	TIS-140	#60 (0.29 mm)	85,700
S21-643-721	TIS-145	#70 (0.25 mm)	98,300
S21-643-725	TIS-150	#80 (0.20 mm)	100,800
S21-643-729	TIS-155	#100 (0.154 mm)	108,400
S21-643-733	TIS-160	#120 (0.132mm)	133,600
S21-643-737	TIS-165	#150 (0.110mm)	136,100
S21-643-741	TIS-170	#180 (91 um)	181,400
S21-643-745	TIS-175	#200 (77 um)	181,400
S21-643-749	TIS-180	#250 (62 um)	272,400
S21-643-753	TIS-185	#300 (45 um)	340,400
S21-643-757	TIS-190	#350 (43 um)	542,200
S21-643-761	TIS-195	#400 (34 um)	605,200

• Size : 1 × 1 m

• Tokyo Screen