

Ro-Tap® Test Sieve Shaker 로텡 체진동기

Cat. No.	Model	적용 Sieve	EA Price ₩
S19-55-250	RX-29-16	Ro-Tap 체진동기 (ϕ 200, ϕ 203mm 용)	7,411,700
S19-55-255	RX-30-16	Ro-Tap 체진동기 (ϕ 300, ϕ 304mm 용)	8,479,300

- Two Dimension Operation의 체진동기
- 수직과 원형운동이 동시에 작동
- Digital Timer 장착 : 99.0 min
- Tyler

Cat. No.	S19-55-250	S19-55-255
Model	RX-29-16	RX-30-16
Oscillation Per Minute	278 \pm 10	
Oscillation Displacement	1-1/8" \times 7/16"	
Taps Per Minute	150 \pm 10	
Testing Application	2" down through 635 mesh (20 μ m)	4" Down Through 500 mesh (25 μ m)
적용 Sieve	ϕ 203 또는 200 mm 최대 6개 (H50 mm)	ϕ 300 또는 304 mm 최대 4개 (H50 mm)
Size	w710 \times d530 \times h40 mm, 82 kg	

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
S19-55-925	10121	Maintenance Kit	305,500
S19-55-935	10122	Repair Kit for only RX-29	1,422,300
S19-55-945	R-30050KR	Sound Enclosure	2,523,900
S19-55-955	H-40041	Test Stand	1,004,600
S19-55-965	LA11112	Wet Sieving Kit for only RX-29	2,865,700

- Maintenance Kit (10121) : 2 Eccentric Discs, 1 Flange Bushing, 1 Bearing, 1 Cork and 2 Flange Bearings
- Repair Kit (10122) : 1 Maintenance Kit, 1 Timing Belt, 1 Upper Carrying Plate, 1 Cam Gear, and 1 Timing Belt Pulley
- R-30050KR : 양카볼트 4개 포함



RX-29-16



Sound Enclosure

Ro-Tap® Test Sieve Shaker 로텡 체진동기

Cat. No.	Model	Description	EA Price ₩
S19-55-270	RX-94-1	Ro-Tap 체진동기 2열	12,272,600

- RX-29, RX-30의 디자인을 동일하게 적용한 체진동기
- 2개의 기둥에 고정된 8" Diameter Test Sieve를 동시에 진동시킴
- Testing Application: 2" Down Through 635 Mesh (20 Micron).
- Clamping System으로 단단히 고정시킬 수 있음
- Tyler

Cat. No.	S19-55-270
Model	RX-94-1
Oscillation Per Minute	278 \pm 10
Oscillation Displacement	1-1/8" \times 7/16"
Taps Per Minute	150 \pm 10
Hammer Weight	5 lbs., 8 oz.
Testing Applications	2" down through 635 mesh (20 micron)
최대 거치 수량	ϕ 203 또는 200mm, 최대 12개 (1열 6개 \times 2열)
Size	w750 \times d680 \times h560 mm, 110 kg

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
S19-55-925	10121	Maintenance Kit	305,500
S19-55-927	10123	Repair Kit for RX-94	1,173,100
S19-55-947	RX-99	Sound Encloser for RX-94	3,604,800
S19-55-967	LA11194	Wet Sieving Kit	3,088,000



- Maintenance Kit (10121)
- 2 Eccentric Discs, 1 Flange Bushing, 1 Bearing, 1 Cork and 2 Flange Bearings
- Repair Kit (10123)
- 1 Maintenance Kit, 1 Timing Belt, 1 Cam Gear, and 1 Timing Belt Pulley