

## I Tall Form Berzelius Beaker, Kimble®

### 장형 유리 비이커

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 적정(Titration) 실험에 알맞음
- 잘 지워지지 않는 Marking Area와 이중 눈금
- 우수한 기계적 강도와 내구성을 가져 화학 약품과 열충격에 강함
- DWK Lifescience / Kimble®

#### without Spout ①

Cat. No.	Model	용량	Size ( $\varnothing \times h, mm$ )	EA Price ₩	PK Price ₩
B05-137-170	14020-100	100 ml	50×78	4,900	12ea 54,800
B05-137-175	14020-200	200 ml	57×104	5,800	12ea 64,500
B05-137-180	14020-300	300 ml	64×118	9,400	12ea 104,500
B05-137-185	14020-400	400 ml	70×127	8,500	6ea 47,300
B05-137-190	14020-600	600 ml	80×150	8,400	6ea 46,700
B05-137-195	14020-1000	1 L	90×188	19,200	6ea 107,600

- Design Meets ASTM Specification E960, Type IV

#### with Spout ②

Cat. No.	Model	용량	Size ( $\varnothing \times h, mm$ )	EA Price ₩	PK Price ₩
B05-137-200	14030-100	100 ml	50×78	5,900	12ea 65,900
B05-137-205	14030-200	200 ml	57×104	6,500	12ea 73,000
B05-137-210	14030-300	300 ml	64×118	8,200	12ea 92,000
B05-137-215	14030-400	400 ml	70×127	9,600	6ea 53,600
B05-137-220	14030-600	600 ml	80×150	13,500	6ea 75,300
B05-137-225	14030-1000	1 L	90×188	19,800	6ea 111,000

- Design Meets ASTM Specification E960, Type III



## I Beaker-Flask Combination, Kimble®

### 플라스크형 비이커, BEAKERplus™

Cat. No.	Model	용량	Size ( $\varnothing \times h, mm$ )	EA Price ₩	PK Price ₩
B05-137-230	14035-150	150 ml	70×83	62,900	6ea 352,700
B05-137-235	14035-250	250 ml	70×107	67,400	6ea 377,600
B05-137-240	14035-500	500 ml	70×178	74,600	6ea 418,000
B05-137-245	14035-1000	1 L	95×222	90,200	2ea 168,600
B05-137-250	14035-1200	1.2 L	95×254	96,900	2ea 181,000

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Griffin Beaker와 삼각 플라스크 결합형
- 읽기 쉬운 흰색 눈금
- 좁게 들어간 형태의 Neck은 비이커를 흔들 때 내용물이 튀는 것을 최소화함
- 깔때기 또는 필터와 사용하거나 미싱용으로 알맞음
- DWK Lifescience / Kimble®



## I Glass Beaker with Handle, Simax®

### 유리 핸들 비이커

Cat. No.	Model	용량	Size ( $\varnothing \times h, mm$ )	EA Price ₩	PK Price ₩
B05-985-337	417011250	250 ml	70×95	6,600	6ea 36,900
B05-985-341	417011400	400 ml	80×110	7,000	6ea 39,000
B05-985-345	417011600	600 ml	90×125	7,700	6ea 42,800
B05-985-349	417011940	1 L	105×145	9,700	2ea 18,000

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Complies with ISO3819, DIN12331
- KAVALIERGLASS

