


CTFE Spoon CTFE 테프론 스푼

Cat. No.	Model	길이	EA Price ₩
S16-64-151	F5013-01	150 mm	33,100
S16-64-155	F5013-02	180 mm	33,900
S16-64-159	F5013-03	210 mm	37,100
S16-64-163	F5013-04	240 mm	40,200
S16-64-167	F5013-05	300 mm	58,300

- 재질 : CTFE (Chlorotrifluoroethylene)
- Flon Industry



CTFE with Spatula Spoon CTFE 테프론 스파츄라 스푼

Cat. No.	Model	길이	EA Price ₩
S16-64-222	F5014-01	150 mm	33,100
S16-64-226	F5014-02	180 mm	33,900
S16-64-230	F5014-03	210 mm	37,100
S16-64-234	F5014-04	240 mm	40,200
S16-64-238	F5014-05	300 mm	58,300

- 재질 : CTFE (Chlorotrifluoroethylene)
- Flon Industry



Disposable Spatula Spoon 일회용 스파츄라 스푼

Cat. No.	Model	Blade	전체 길이	Description	PK Price ₩
S16-100-199	6284543	w17 mm	178 mm	Non-sterilized	100ea 49,000
S16-100-201	6284544			Sterilized	100ea 122,300

- 재질 : PS, 내화학성, 내부식성
- Clean Room Class 7의 환경에서 생산 및 포장
- Complies with EU Food requirements, FDA regulation
- Wenk Labtec



Plastic Spoon 플라스틱 스푼

- 재질 : Impact Resistant PS
- Brand GmBH

Double Spoon Type ①

Cat. No.	Model	길이	End Size	EA Price ₩
S16-59-105	139810	150 mm	30 mm	4,100
S16-59-110	139815	180 mm		4,200

Spatula / Spoon Type ②

Cat. No.	Model	길이	End Size	EA Price ₩
S01-59-115	139817	180 mm	30 mm	4,300
S01-59-120	139820	210 mm		4,700



Porcelain Spatula with Scoop 자제 스파츄라

Cat. No.	Model	길이	EA Price ₩
S16-128-103	218/1	103 mm	6,600
S16-128-105	218/1a	148 mm	7,000
S16-128-107	218/1b	164 mm	7,800
S16-128-109	218/2	187 mm	9,800
S16-128-111	218/2a	200 mm	12,300
S16-128-113	218/3	250 mm	13,600

- Max. Temp. : 1000°C
- Jipo / LabCeramic®

