

Aluminum Scoop 알루미늄 스코프

Cat. No.	Model	용량	전체 길이	EA Price ₩
S16-102-511	3720	60 ml	150 mm	5,900
S16-102-518	3721	130 ml	200 mm	9,900
S16-102-525	3722	250 ml	250 mm	12,000
S16-102-532	3724	450 ml	300 mm	19,500
S16-102-539	3725	800 ml	330 mm	33,300
S16-102-546	3726	1.5 L	390 mm	50,200
S16-102-553	3727	2.5 L	410 mm	65,200

- 재질 : Aluminum
- LK Labkorea



Polycarbonate Scoop 투명 PC 스코프

Cat. No.	Model	용량	Size	EA Price ₩
S16-517-015	PCS015	150 ml	w60×L242 mm	6,000
S16-517-030	PCS030	250 ml	w80×L267 mm	6,000
S16-517-060	PCS060	500 ml	w105×L316 mm	6,900
S16-517-090	PCS090	1 L	w126×L334 mm	9,000
S16-517-150	PCS150	1.5 L	w157×L358 mm	10,900

- 재질 : PC (투명), 충격 및 열에 강한 스코프
- 사용온도 : -38 ~ 100°C
- NSF 인증
- Rayhope®



Slope-type Scoop, PE 경사형 스코프

Cat.No.	Model	용량	전체 길이	EA Price ₩
S16-59-321	139035	500 ml	350 mm	10,600
S16-59-325	139040	1 L	400 mm	15,100

- 재질 : PE, -10 ~ 110°C 내열
- Brand GmBH



PE Bio Clean Scoop 바이오 친환경 스코프

- 친환경 PE로 제조되어 재활용 가능
- Clean Room Class 7의 환경에서 생산 및 개별 포장
- Complies with EU Food requirements, FDA regulation
- Wenk Labtec

Non-sterilized

Cat. No.	Model	용량	전체 길이	PK Price ₩
S16-100-131	6265055	25 ml	141 mm	10ea 43,500
S16-100-134	6265056	50 ml	160 mm	10ea 44,900
S16-100-137	6265057	100 ml	205 mm	10ea 46,400
S16-100-140	6265058	150 ml	216 mm	10ea 47,600

Sterilized

Cat. No.	Model	용량	전체 길이	PK Price ₩
S16-100-143	6265049	25 ml	141 mm	10ea 48,400
S16-100-146	6265050	50 ml	160 mm	10ea 49,900
S16-100-149	6265051	100 ml	205 mm	10ea 51,400
S16-100-152	6265052	150 ml	216 mm	10ea 52,800

