

Plastic Scoop, PP 플라스틱 스코프

Cat. No.	Model	용량	Scoop	전체 길이	EA Price ₩
S16-20-121	139102	2 ml	w12×L33	60 mm	700
S16-20-125	139103	5 ml	w15×L45	80 mm	800
S16-20-129	139105	10 ml	w18×L53	100 mm	1,000
S16-20-133	139110	25 ml	w22×L73	135 mm	1,200
S16-20-137	139115	50 ml	w30×L90	160 mm	2,200
S16-20-141	139120	100 ml	w37×L115	200 mm	3,300
S16-20-145	139125	250 ml	w50×L155	260 mm	6,600
S16-20-149	139130	500 ml	w60×L195	315 mm	9,300
S16-20-153	139135	1 L	w80×L247	385 mm	11,000

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Self-Standing, Flat Top
- 고체, 액체 모두 사용 가능하며 눈금 표시
- Brand GmBH



Polypropylene Scoop 플라스틱 스코프, Antistatic PP

Cat. No.	Model	Size	Weight	EA Price ₩
S16-100-161	6901785	w110×d150×L250 mm	90 g	16,100
S16-100-163	7633782	w138×d187×L310 mm	187 g	24,900
S16-100-165	7618436	w167×d230×L355 mm	289 g	29,500


- 재질 : PP, -10 ~ 135°C 내열, 내화학성, Autoclavable
- 제약, 생명공학, 식품, 사료, 화학, 산업, 농업, 양식업, 샘플링용으로 사용 가능
- 세척이 쉽고 매끄러우며 Pore-free 표면으로 식품 분야에 적합
- Complies with EU Food requirements, FDA regulation
- Wenk Labtec



PTFE Scoop PTFE 스코프

Cat. No.	Model	전체 길이	Size	EA Price ₩
S16-450-333	041.001	140 mm	w30×h80×Ø20 mm	44,300

- 재질 : PTFE, -200 ~ 280°C 내열
- 불투명 PTFE 테프론 재질로 최상의 내화학성 가짐
- Cowie



Stainless Steel Scoop 스테인레스 스코프

Cat. No.	Model	용량	Scoop (Ø×h)	전체 길이	EA Price ₩
S16-459-013	SCP050	500 ml	98×150 mm	270 mm	27,700
S16-459-016	SCP100	1 L	115×160 mm	280 mm	31,700
S16-459-019	SCP200	2 L	130×180 mm	300 mm	35,700
S16-459-023	SCP300	3 L	150×220 mm	340 mm	39,600

- 재질 : Stainless Steel 304
- Made in Korea
- LK Labkorea



Deflagration Spoon 연소 스푼

Cat. No.	Model	길이	EA Price ₩
S16-139-841	1270-20	200 mm	3,300
S16-139-843	1270-30	300 mm	4,500
S16-139-845	1270-40	400 mm	5,400

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- LK Labkorea

