

Scoop / Ladle 국자

• Usbeck

Sodium Spoon, Nickel Plated Brass

| Cat. No. | Model | Scoop | 전체 길이 | EA Price ₩ |
|-------------|-------|--------|---------|------------|
| S16-284-101 | 3411 | Ø25 mm | L180 mm | 13,000 |

Deflagrating Spoon, Stainless Steel 18/10

| Cat. No. | Model | Scoop | 전체 길이 | EA Price ₩ |
|-------------|-------|--------|---------|------------|
| S16-284-105 | 3420 | Ø16 mm | L450 mm | 5,200 |

Ladle, Stainless Steel 18/10

| Cat. No. | Model | Scoop | 전체 길이 | EA Price ₩ |
|-------------|-------|--------|---------|------------|
| S16-284-109 | 3430 | Ø75 mm | L275 mm | 9,400 |



PFA Small Size Dipper PFA 테프론 국자

| Cat. No. | Model | 용량 | Dipper (Ø x h) | 전체 길이 | EA Price ₩ |
|-------------|----------|--------|----------------|--------|------------|
| S16-102-920 | F1212-01 | 180 ml | 86 x 45 mm | 400 mm | 501,900 |


• 재질 : PFA • Flon Industry



PFA Dipper & Sieve PFA 테프론 망 국자 ①

| Cat. No. | Model | 용량 | Dipper (Ø x h) | Mesh | EA Price ₩ |
|-------------|----------|-------|----------------|--------|------------|
| S16-102-956 | F1104-02 | 50 ml | 69 x h22 mm | 1.0 mm | 159,800 |
| S16-102-959 | F1104-03 | | | X | 157,300 |
| S16-460-250 | F997-03 | 80 ml | 75 x h19 mm | 1.0 mm | 214,500 |

• 재질 : PFA
• F1104-02, F1104-03 전체 길이 : L175 mm
• Flon Industry / BOLA (F997-03)



PTFE Dipper Vessel PTFE 테프론 디퍼 용기 ②

| Cat. No. | Model | 용량 | 외부 Size | I.D. | EA Price ₩ |
|-------------|----------|--------|-----------------|--------|------------|
| W07-460-011 | H1138-08 | 90 ml | Ø 38 x h 100mm | Ø 35mm | 302,000 |
| W07-460-015 | H1138-16 | 240 ml | Ø 60 x h 175mm | Ø 57mm | 370,800 |
| W07-460-019 | H1138-24 | 500 ml | Ø 100 x h 230mm | Ø 95mm | 608,900 |

• 재질 : PTFE, -200 °C to 250°C 내열 • Wall과 Bottom에 Hole이 있는 용기
• 용도 : Aggressive or Pure Substance안에 시료를 Dipping & Washing & Rinsing
• BOLA

Ladle Scoop 국자

• 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
• LK Labkorea

일자형, Flat Handle

| Cat. No. | Model | 용량 | Scoop | 전체 길이 | EA Price ₩ |
|-------------|-------|--------|---------|--------|------------|
| S16-102-560 | 8740 | 90 ml | Ø70 mm | 185 mm | 33,100 |
| S16-102-567 | 8741 | 200 ml | Ø90 mm | 315 mm | 39,900 |
| S16-102-574 | 8742 | 220 ml | Ø100 mm | | 46,500 |

곡형, Round Handle

| Cat. No. | Model | 용량 | Scoop | 전체 길이 | EA Price ₩ |
|-------------|-------|--------|---------|--------|------------|
| S16-102-581 | 8791 | 100 ml | Ø70 mm | 300 mm | 66,500 |
| S16-102-588 | 8792 | 175 ml | Ø80 mm | | 65,400 |
| S16-102-595 | 8794 | 250 ml | Ø100 mm | 400 mm | 91,400 |

