

## Micro Absorption Cell, 2-Side Polished 마이크로 흡광 셀, 2면 투명

- 재질 : Spectrosil® Far UV Quartz
- 투과 범위 : 170 ~ 2700 nm
- PTFE 재질의 Cover 또는 Stopper 포함
- 표준 셀 또는 세미 마이크로 셀 보다 적은 양의 샘플을 사용할 때 적합
- Starna

### Type 18 (Clear Wall), with PTFE Cover

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-132	18-Q-5	5 mm	0.35 ml	12.5×7.5×45	2 mm	273,900
S18-77-136	18-Q-10	10 mm	0.7 ml	12.5×12.5×45		212,200
S18-77-140	18-Q-20	20 mm	1.4 ml	12.5×22.5×45		339,300
S18-77-144	18-Q-40	40 mm	2.8 ml	12.5×42.5×45		427,000
S18-77-148	18-Q-50	50 mm	3.5 ml	12.5×52.5×45		471,100
S18-77-152	18-Q-100	100 mm	7.0 ml	12.5×102.5×45		747,300
S18-77-156	18/9-Q-10	10 mm	0.58 ml	12.5×12.5×45		265,400

- Base Thickness : 3 mm (단, 18/9-Q-10은 9 mm)

### Type 18B (Black Wall), with PTFE Cover

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-732	18B-Q-5	5 mm	0.35 ml	12.5×7.5×45	2 mm	481,600
S18-77-736	18B-Q-10	10 mm	0.7 ml	12.5×12.5×45		384,200
S18-77-740	18B-Q-20	20 mm	1.4 ml	12.5×22.5×45		576,100
S18-77-744	18B-Q-40	40 mm	2.8 ml	12.5×42.5×45		726,300
S18-77-748	18B-Q-50	50 mm	3.5 ml	12.5×52.5×45		863,400
S18-77-756	18B/9-Q-10	10 mm	0.58 ml	12.5×12.5×45		498,700

- Base Thickness : 3 mm (단, 18B/9-Q-10은 9 mm)

### Type 28 (Clear Wall), with PTFE Stopper

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-160	28-Q-5	5 mm	0.35 ml	12.5×7.5×48	2 mm	375,300
S18-77-164	28-Q-10	10 mm	0.7 ml	12.5×12.5×48		310,700
S18-77-168	28-Q-20	20 mm	1.4 ml	12.5×22.5×48		466,100
S18-77-172	28-Q-40	40 mm	2.8 ml	12.5×42.5×48		585,500
S18-77-176	28-Q-50	50 mm	3.5 ml	12.5×52.5×48		645,700
S18-77-180	28/9-Q-10	10 mm	0.58 ml	12.5×12.5×48		353,900

- Base Thickness : 3 mm (단, 28/9-Q-10은 9 mm)

- Light Path 40, 50 mm 제품은 Stopper 2개 포함

### Type 28B (Black Wall), with PTFE Stopper

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-760	28B-Q-5	5 mm	0.35 ml	12.5×7.5×48	2 mm	661,800
S18-77-764	28B-Q-10	10 mm	0.7 ml	12.5×12.5×48		546,100
S18-77-768	28B-Q-20	20 mm	1.4 ml	12.5×22.5×48		819,000
S18-77-772	28B-Q-40	40 mm	2.8 ml	12.5×42.5×48		1,032,000
S18-77-776	28B-Q-50	50 mm	3.5 ml	12.5×52.5×48		1,136,000
S18-77-780	28B/9-Q-10	10 mm	0.58 ml	12.5×12.5×48		619,000

- Base Thickness : 3 mm (단, 28B/9-Q-10은 9 mm)

- Light Path 40, 50 mm 제품은 Stopper 2개 포함

