

Anaerobic Absorption Cell with Screw Cap, 2-Side Polished

스크류 캡 흡광 셀, 2면 투명

- 재질 : Spectrosil® Far UV Quartz
- 투과 범위 : 170 ~ 2700 nm
- GL14 Septum Cap (with 9 mm Hole) 또는 Closed Cap 포함
- Cap Max. Pressure : 5×10^5 pa (5 bar)
- 혐기조건 하에서 사용 가능
- Starna

with Septum Cap

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-256	1GL14-Q-1-S	1 mm	0.4 ml	12.5×3.5×66	10 mm	437,000
S18-77-260	1GL14-Q-2-S	2 mm	0.8 ml	12.5×4.5×66		437,000
S18-77-264	1GL14-Q-5-S	5 mm	1.6 ml	12.5×7.5×66		421,100
S18-77-268	1GL14-Q-10-S	10 mm	3.5 ml	12.5×12.5×66	4 mm	375,300
S18-77-272	9GL14-Q-10-S		1.4 ml			437,000
S18-77-276	9BGL14-Q-10-S		0.7 ml			571,200
S18-77-280	18GL14-Q-10-S	10 mm	0.7 ml	12.5×12.5×66	2 mm	447,500
S18-77-284	18BGL14-Q-10-S					611,600

- 9BGL14-Q-10-S & 18BGL14-Q-10-S : Black Cell


with Closed Cap

Cat. No.	Model	Light Path	용량	외부 Size (w×d×h, mm)	내부 Width	EA Price ₩
S18-77-258	1GL14-Q-1-C	1 mm	0.4 ml	12.5×3.5×66	10 mm	437,000
S18-77-262	1GL14-Q-2-C	2 mm	0.8 ml	12.5×4.5×66		437,000
S18-77-266	1GL14-Q-5-C	5 mm	1.6 ml	12.5×7.5×66		421,100
S18-77-270	1GL14-Q-10-C	10 mm	3.5 ml	12.5×12.5×66	4 mm	375,300
S18-77-274	9GL14-Q-10-C		1.4 ml			437,000
S18-77-278	9BGL14-Q-10-C		0.7 ml			571,200
S18-77-282	18GL14-Q-10-C	10 mm	0.7 ml	12.5×12.5×66	2 mm	447,500
S18-77-286	18BGL14-Q-10-C					611,600

- 9BGL14-Q-10-C & 18BGL14-Q-10-C : Black Cell



9BGL14-Q-10



Short Light Path Absorption Cell, 2-Side Polished

단파 흡광 셀, 2면 투명

- 재질 : Spectrosil® Far UV Quartz
- 투과 범위 : 170 ~ 2700 nm
- 분해 가능
- 파장이 단파장이거나 크기가 큰 Window 표면이 필요할 때 사용되는 제품
- 불투명 액체 또는 비스코스(Viscose) 액체 사용시 유용
- CH/2049 Cell Holder와 함께 사용
- Starna

Type 20/O, Open-Ended

Cat. No.	Model	Light Path	용량 (ml)	외부 Size (w×d×h, mm)	내부 Size (w×h)	EA Price ₩
S18-77-288	20/O-Q-0.01	0.01 mm	0.004	12.5×2.5×45	10×43.5 mm	420,800
S18-77-292	20/O-Q-0.05	0.05 mm	0.020	12.5×2.5×45		401,600
S18-77-296	20/O-Q-0.1	0.1 mm	0.040	12.5×2.6×45		336,700
S18-77-300	20/O-Q-0.2	0.2 mm	0.080	12.5×2.7×45		336,700
S18-77-304	20/O-Q-0.5	0.5 mm	0.190	12.5×3.0×45		326,300

Type 20/C, Closed

Cat. No.	Model	Light Path	용량 (ml)	외부 Size (w×d×h, mm)	내부 Size (w×h)	EA Price ₩
S18-77-308	20/C-Q-0.008	0.008 mm	0.002	12.5×2.508×45	8×38 mm	1,364,800
S18-77-312	20/C-Q-0.01	0.01 mm	0.003	12.5×2.51×45		420,800
S18-77-316	20/C-Q-0.05	0.05 mm	0.015	12.5×2.55×45		401,600
S18-77-320	20/C-Q-0.1	0.1 mm	0.030	12.5×2.6×45		347,300
S18-77-324	20/C-Q-0.2	0.2 mm	0.060	12.5×2.7×45		347,300
S18-77-328	20/C-Q-0.5	0.5 mm	0.150	12.5×3.0×45		347,300

Cell Holder

Cat. No.	Model	적용 셀	EA Price ₩
S18-203-A20	CH/2049	Light Path 5 mm 이하 사각 Cell용	131,800

