



500 ml Centrifuge Bottle			
Sl	Component	Material	Qty
1	Container	PP	1
2	Closure	PP	1

CUSTOMER COPY
DATE..... 04.05.2023

PP/PP	---	96±10 gms	Ø69.3X166mm	450 ml	00	---	---	---
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name:			
DRAWN	Arunopal		03.05.2023		Centrifuge Bottle 500 ml			
CHECKED	Iranna		03.05.2023					
APPROVED	Suresh		03.05.2023					

Drawing Number.		A4
544020		
SCALE:NTS		R00
SHEET 1 OF 1		

Tarsons
TRUST DELIVERED

Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Products Limited.

Unspecified tolerances for open dimensions	
0.00	± 0.05 mm
0.0	± 0.1 mm
0	± 0.5 mm
Angular	± 0.2°