



1000 ml Centrifuge Bottle		
Sl	Component	Material
1	Container	PP
2	Closure	PP

**CUSTOMER COPY**  
DATE..... 05.05.2023

PP/PP	---	155±20 gms	Ø96.7X167mm	915 ml	00
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name: <b>Centrifuge Bottle 1000 ml (Short Neck)</b>
DRAWN	Arunopal		03.05.2023		
CHECKED	Iranna		03.05.2023		
APPROVED	Suresh		03.05.2023		

Revised Details	Revised by	Date
Drawing Number. <b>544040</b>		
SCALE:NTS		SHEET 1 OF 1

**Tarsons**  
TRUST DELIVERED

Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Products Limited.

Unspecified tolerances for open dimensions	
0.00	± 0.05 mm
0.0	± 0.1 mm
0	± 0.5 mm
Angular	± 0.2°

A4  
R00