



SL No	Description	Material	QTY
1	Bottle	HDPE	1
2	Closure	PP	1

CUSTOMER COPY
DATE: 23-08-2022

Material	Shrinkage	Approx Weight	Ø93.1x214	Brimfull volume approx	00	customer copy drawing	Arunopal	23.08.2022
	NAME	SIGNATURE	DATE	Angle of Projection	Rev.no	Revised Details	Revised by	Date
DRAWN	Arunopal		23.08.2022		Part Name: 1000 ml Round Narrow Mouth Bottle(Amber)			
CHECKED	Iranna							
APPROVED	Suresh							

Tarsons Products Ltd.

Remarks:
This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Product Limited.

Unspecified tolerances for open dimensions	
0.00	± 0.05 mm
0.0	± 0.1 mm
0	± 0.5 mm
Angular	± 0.2°

Drawing Number.	A4
581250	R00
SCALE:1:1	SHEET 1 OF 1