


If in doubt please ask
All dimensions are in mm



CUSTOMER COPY
DATE 08-04-2022

PP	%	36.7 ± 5.0 gms		wxbxh		Sl no	Description	Revisions	Date
MATERIAL	Shrinkage	WEIGHT		FINISH SIZE	Brimfull volume :	DO NOT SCALE DRAWING		Revision	
	NAME	SIGNATURE	DATE	TITLE:	Angle of Projection	Description 250 ml Round Wide Mouth Bottle			
DRAWN	Shubham		08-04-2022						
CHK'D	Iranna				Remarks:				
APPVD	Suresh								
DEBURR AND BREAK SHARP EDGES		Unspecified tolerances for open dimension				DWG NO.		A4	
		0 - 50mm	+/-0.1 mm	This is confidential document. It must not be copied wholly or in part not be disclosed to third party or not be amended without prior written consent from M/S tarsons products pvt ltd		582230			
		50 - 100mm	+/-0.2mm						
		100 <	+/-0.25mm						
		Angular	+/- 0.5°			SCALE:NTS		SHEET 1 OF 1	

