

《 규 격 서 》

카달로그 번호	품목번호	품목 및 규격	단위	수량
Cat. No	Item No	Description	Unit	Quantity
B04-01-050	COB10	Circulating Oil Bath	set	1

1 Feature

- 내부 순환기능을 이용하여 Bath 내부 온도를 균일하게 유지
- 강력한 순환 펌프를 이용하여 외부기기와 연결하여 사용하기에 알맞음
- 10단계 펌프 유량 조절
- 0.01°C 정밀 제어 구현
- 요청시 Cooling Coil 장착 가능
고온 구동 후 외부에서 냉각수를 순환시켜 빠르게 온도를 떨어트릴 수 있음 (베이스 커버에 장착)
- Bath 내부 용액을 편리하게 배출 시킬 수 있는 Drain 장착
- 외부 센서 및 RS-485 PC 통신 기능 제공 (통신프로그램은 엘케이랩코리아 홈페이지에서 다운로드)
- Graphic LCD Display, 간편하고 직관적인 Controller 조작
- Auto-tuning 기능 및 온도 Calibration 기능 내장
- 2 가지 타입의 타이머. 구동과 동시에 시작, 목표 온도 도달 후 시작
- Automatical Fault Sensing System 기기 스스로 이상 상황을 감시 (특허 등록)
- Low Water Alarm, Bath 용액 부족으로 인한 Heater 과열방지
- 2중화 과온 안전장치, 1차 Controller 과온 Alarm과 2차 과온 차단 회로

2 Specifications

Cat. No.		B04-01-050
Model		COB10
Capacity		10 L
Controller	Control	PID Control, Auto-tuning
	Display	GLCD (Graphic LCD)
	Resolution	0.01°C
	Timer	Run Start / Wait start (99 hr 59 min)
Temperature	Range	주변온도 +10°C to +250°C
	Accuracy at 100°C	±0.3°C
	Uniformity at 100°C	±0.42°C
Pump	Capacity	22 L/min (at water 25°C)
	Pressure	0.47 bar / 6.8 psi (at water 25°C)

Size	Internal	w240×d300×h150 mm
	External	w310×d370×h425 mm
	Top Open	w240×d130 mm
Material	Interior	Stainless Steel 304
	Exterior	Powder Coated Steel
Electric Supply	Power	1 Phase, 220 VAC, 50/60Hz
	Max Consumption	1.1 kW (5 A)
	Power Line	Standard Plug
Other	Drain	1/4" Hose Nipple Male
	Circulating IN / OUT	1/2" Stainless Steel Tubing
Safety		Over Temp. & Current Protection, Automatic Fault Sensing System

3 Accessory

- B04-01-085 / LI-USB / RS485 to USB Communication Cable
- B04-91-483 / ET101M / External Sensor, Pt100, Ø4xL300mm
- C20-184-035 / S050-U / 1/2" Stainless Steel Tubing Connector

4 Remark

- 사용자의 지정장소에 설치 후 1년 동안 무상보증.
- PL Insurance (제조물배상책임보험)
- Validation Service (IQ, OQ, PQ) 가능