

《 규 격 서 》

| 카달로그 번호 | 품목번호 | 품목 및 규격 | 단위 | 수량 |
|------------|---------|----------------------|------|----------|
| Cat. No | Item No | Description | Unit | Quantity |
| B04-01-070 | COB30 | Circulating Oil Bath | set | 1 |

1 Feature

- 내부 순환기능을 이용하여 Bath 내부 온도를 균일하게 유지
- 강력한 순환 펌프를 이용하여 외부기기와 연결하여 사용하기에 알맞음
- 10단계 펌프 유량 조절
- 0.01°C 정밀 제어 구현
- 요청시 Cooling Coil 장착 가능
고온 구동 후 외부에서 냉각수를 순환시켜 빠르게 온도를 떨어트릴 수 있음 (베이스 커버에 장착)
- Bath 내부 용액을 편리하게 배출 시킬 수 있는 Drain 장착
- 외부 센서 및 RS-485 PC 통신 기능 제공 (통신프로그램은 엘케이랩코리아 홈페이지에서 다운로드)
- Graphic LCD Display, 간편하고 직관적인 Controller 조작
- Auto-tuning 기능 및 온도 Calibration 기능 내장
- 2 가지 타입의 타이머. 구동과 동시에 시작, 목표 온도 도달 후 시작
- Automatic Fault Sensing System 기기 스스로 이상 상황을 감시 (특허 등록)
- Low Water Alarm, Bath 용액 부족으로 인한 Heater 과열방지
- 2중화 과온 안전장치, 1차 Controller 과온 Alarm과 2차 과온 차단 회로

2 Specifications

| | | |
|-------------|---------------------|---------------------------------------|
| Cat. No. | | B04-01-070 |
| Model | | COB30 |
| Capacity | | 30 L |
| Controller | Control | PID Control, Auto-tuning |
| | Display | GLCD (Graphic LCD) |
| | Resolution | 0.01°C |
| | Timer | Run Start / Wait start (99 hr 59 min) |
| Temperature | Range | 주변온도 +10°C to +250°C |
| | Accuracy at 100°C | ±0.3°C |
| | Uniformity at 100°C | ±1.23°C |
| Pump | Capacity | 22 L/min (at water 25°C) |
| | Pressure | 0.47 bar / 6.8 psi (at water 25°C) |

| | | |
|-----------------|----------------------|--|
| Size | Internal | w300×d450×h230 mm |
| | External | w370×d520×h505 mm |
| | Top Open | w320×d270 mm |
| Material | Interior | Stainless Steel 304 |
| | Exterior | Powder Coated Steel |
| Electric Supply | Power | 1 Phase, 220 VAC, 50/60Hz |
| | Max Consumption | 2.8 kW (12.8 A) |
| | Power Line | Standard Plug |
| Other | Drain | 1/4" Hose Nipple Male |
| | Circulating IN / OUT | 1/2" Stainless Steel Tubing |
| Safety | | Over Temp. & Current Protection, Automatic Fault Sensing System |

3 Accessory

- B04-01-085 / LI-USB / RS485 to USB Communication Cable
- B04-91-483 / ET101M / External Sensor, Pt100, Ø4xL300mm
- C20-184-035 / S050-U / 1/2" Stainless Steel Tubing Connector

4 Remark

- 사용자의 지정장소에 설치 후 1년 동안 무상보증.
- PL Insurance (제조물배상책임보험)
- Validation Service (IQ, OQ, PQ) 가능