

# 《 규 격 서 》

카달로그 번호	품목번호	품목 및 규격	단위	수량
Cat. No	Item No	Description	Unit	Quantity
C10-91-203	LM-0900P	Desktop Clean Hood	set	1

## 1 Feature

- 우기미량 분석용 시험장비의 준비작업 등 청정공간을 유지하여 약품 증기 유출을 차단하는 장비
- 청정 공간에서 산 작업이 필요한 경우에 사용
- 공기가 상부에서 필터를 거쳐 수직 방향으로 하강하여 후면으로 배기 되는 방식
- 내부가 비금속 재질로 구성되었기 때문에 부식에 의한 오염이 없음
- 균일한 Lamina Flow 를 형성하여 신뢰성을 높임
- Class 100 HEPA Filter를 사용하여 99.97%의 집진율
- FL Lamp 및 UV Lamp 적용
- Filter 교환시기를 경보음과 LCD 화면에 표시
- 차압계로 필터의 상태를 확인 할 수 있음
- 10단계 풍속 조절 가능
- 2단 접이식 홀딩 도어
- Timer 기능 (99 hour 59 min)

## 2 Specifications

Cat. No.		C10-91-203
Model		LM-0900P
Size	Internal	w880×d570×h590 mm
	External	w900×d720×h1000 mm
	Exhausted Duct	L150, d200 mm
Filter	HEPA	HEPA filter 0.3 um 99.97% Efficiency Filter
	Pre-filter	Efficiency Nylon Filter
Material	Upper Case	Powder Coated Steel
	Working Space	PVC
Wind Velocity		0.3 ~ 0.5 m/sec, 10 step adjustable
Lamp		FL lamp 1ea, UV lamp 1ea
Control System		LCD Digital Controller
Power		1 Phase, 220 Vac, 60 Hz

### 3 Accessory

- C10-91-903 / LM-FL10 / FL lamp 10W
- C10-91-913 / LM-UV10 / UV lamp 10W
- C10-91-941 / LM-CHF900 / Filter for LM-0900P

### 4 Remark

- 사용자의 지정장소에 설치 후 1년 동안 무상보증.
- PL Insurance (제조물배상책임보험)
- Validation Service (IQ, OQ, PQ) 가능