

Data sheet

Product: PORAFIL® TE

REF: 670300090

Membrane material: Polytetrafluoroethylene (PTFE)

Color: white

Air Permeability: 17.5 – 27.5 L/s/m²

Thickness: 0.15 – 0.25 mm

Water Entry Pressure: ≥ 0,15 bar

Temperature stability: max. 145 °C

Pore diameter

0.20 µm 0.45 µm 1.0 µm 3.0 µm

x

Filter diameter

13 mm 25 mm 47 mm 50 mm 90 mm 100 mm 142 mm 220 mm 293 mm

x

These data are typical values and must not be regarded as specification.

