

COMMERCIAL CENTER&R-D Center
 #203, Building 2, 89 Shengli Rd, Jiangning District, Nanjing 211106, China
 Tel.: 0086 25 86216803
 Fax: 0086 25 86214463
 E-mail: info@citotest.com

PRODUCTION
 #339, Beihai Road, Haimen City, Jiangsu 226100, P.R.China
 Tel.: 0086 513 82259348
 Fax: 0086 513 82110195
 E-mail: info@citotest.com.cn



CITOTEST LABWARE MANUFACTURING CO.,LTD

www.citotest.com

PRODUCT DATA SHEET

GENERAL INFORMATION

PRODUCT NAME: Microscope Cover Glass

FURTHER DESCRIPTION: Microscope Cover Glass, Circular, Dia.14mm. Thickness 1.0# 0.13-0.16mm. Economical Package.100pcs/box,10x100pcs/inner box,wrapped by aluminum vacumm pack,500box/carton

P/N: 0344-1420

CATEGORY NUMBER: 10210014CE

SPECIFICATION

COLOR: Transparent

MATERIAL: Super white glass

COMPONENT: SiO₂(73%):Na₂O(14%):K₂O(1.5%):CaO(5.82%):MgO(4%):Al₂O₃(1.5%):SO₃(0.15):Fe₂O₃(0.03%)


DIMENSION (mm): ø14

THICKNESS(mm): 0.13-0.16

PRODUCT PICTURE:



TECHNICAL PROCESS:



INTENDED APPLICATION

To keep solid specimens pressed flat or liquid samples shaped into a flat layer of even thickness on microscope slides.

PACKING

100pcs/box,10x100pcs/inner box,wrapped by aluminum vacumm pack,500box/carton

BOX:



样品编号: 0344-1420 客户编号: 物料代码: 4.07.10444 标签尺寸: 23x23mm
 颜色: PANTONE 534C PANTONE 187C

INNERBOX: [Innerbox Label Size(mm)]



CARTON: [Carton Label Size(mm) 100×130]



STERILIZATION

STERILE EO STERILE R STERILE A 

STORAGE INSTRUCTION

Store in a cool, dry place.

TERM OF VALIDITY

18 months

REMARKS

Conforms with the essential requirements of the European in Vitro Diagnostic 98/79 EC.