

Schlenk O-ring Type Vacuum Trap 오링형 쉬링크 진공 트랩

Cat. No.	Model	Body (od×h, mm)	O-ring	EA Price ₩
V03-132-901	AF-0200-01	57×300	15	476,900

- Borosilicate Glass 3.3
- O-ring Joint형으로 Stopper에 Trap Body 연결
- Two-piece Vacuum Tight, Grease Free, 57mm Flat O-ring Flange
- 구성 : Trap Stopper and Body, #331 Viton o-ring 1개,
#116 Viton o-ring 2개, Horseshoe Clamp
- Chemglass

Replacement Part

Cat. No.	Model	Description	EA Price ₩
V08-132-062	AF-0200-02	Trap Body, #15 O-ring	111,300
V08-132-963	AF-0200-03	Stopper, #15 O-ring	263,700
V08-132-916	CG-139-02	57 mm Horseshoe Clamp	81,000
V08-132-644	CG-305-116	#116 Viton O-ring	2,500
V08-132-914	CG-305-331	#331 Viton O-ring	8,000



Vacuum

O-Ring Type Trap / Bubbler, Pressure Release / Check Valve